INOVOFILL®

Instructions for Use

EN

INOVOFILL®

High molecular weight cross-linked sodium hyaluronate.

DESCRIPTION

- **INOVOFILL®** is a non-pyrogenic, viscoelastic solution consisting of a purified sodium hyaluronate for intra-articular injections ready to use. Sodium hyaluronate in **INOVOFILL®** is manufactured by fermentation and is chemically crosslinked.
- **INOVOFILL®** is produced according to the body's physiological pH and osmolality with no animal protein.
- The administration of **INOVOFILL®** results in an improvement of articular mobility and pain reduction.
- Each pre-filled syringe with 3 mL of **INOVOFILL®** contains:
 - Cross-linked Sodium Hyaluronate: 30 mg
 - Phosphate Buffer pH 7 q.s.; 3 ml

INDICATION

- **INOVOFILL®** is indicated for patients with osteoarthritis (OA) who have failed to respond adequately to conventional non-pharmacological treatments and conventional anti-inflammatories (e.g., ibuprofen and celecoxib).

CONTRAINDICATIONS

- Do not administer to patients with known hypersensitivity (allergy) to hyaluronan (sodium hyaluronate) and its derivatives.
- Do not administer to patients with inflamed or infected in articular fluid.
- **INOVOFILL®** should not be used in patients with skin diseases or infections in the area of the injection or joint.
- The safety and effectiveness of **INOVOFILL®** have not been established in pediatric, pregnant women, and lactating women.
- **INOVOFILL®** has not been tested to show pain relief for conditions other than OA.

PRECAUTIONS

- Only medical professionals trained in the accepted injection techniques for delivering agents into the joint should inject **INOVOFILL®**.
- Not all patients respond equally well to this medical device.
- It is recommended that patients avoid intense or prolonged weight-bearing activities within 48 hours following the injection.
- Patients should be instructed to contact their treating physician if signs of skin disorder (such as change of color or open sores) appear.
- Strict aseptic administration technique must be followed.
- Do not use **INOVOFILL®** if the package has been opened or damaged.
- The pre-filled syringe is intended for single use only. The contents of the syringe should be used immediately after opening, and discard any unused **INOVOFILL®**.
- Remove joint effusion, if present, before injecting **INOVOFILL®**.
- **INOVOFILL®** should be used with caution when there is signs of lymphatic or venous stasis at the injection site.

WARNINGS

- Do not concomitantly use disinfectants containing quaternary ammonium salts (such as benzalkonium chloride) for skin preparation because hyaluronan can precipitate in their presence.
- Do not inject **INOVOFILL®** extra-articular, into the synovial tissues or joint capsule.
- Do not inject **INOVOFILL®** intravascularly, it may cause systemic adverse events.
- The mixing or simultaneous administration with other intra-articular injectables is not recommended.
- Do not use the same syringe for removing synovial fluid and injecting **INOVOFILL®**.

ADVERSE EFFECTS


- Mild to moderate pain, swelling or stiffness may occur after the intra-articular injection. These reactions are generally transient and usually resolved spontaneously or with conservative treatment.

ASSEMBLY OF NEEDLE TO SYRINGE

- For optimal use of **INOVOFILL®**, it is important that the needle is properly connected to the syringe. See diagrams 1, 2, 3 and 4.
- Fig 1, Firmly hold the glass cylinder of the syringe and the Luer-lock adaptor between the thumb and the index fingers.
- Fig 2, Grasp the protective cap with the other hand and break the seal.
- Fig 3, Push & Twist the needle on the syringe until a resistance is felt. Do not over-tight. Overtightening of the needle may lead to the Luer-lock moving and dislodging from the syringe.
- Fig 4, Keep holding the Luer-lock moving and remove the sheath from the needle.

STORAGE

- Store between 2 and 25°C. Protected from light and freezing. Avoid mechanical shocks.

Marketing Authorization Holder:
MEDISA ARA GOSTAR IRANIAN
Number 21, Jooybar Alley Shad St,
Molae Sadra St, Tehran, Iran
Manufacturer:
ARYA TAJHIZ IBN SINIA
Number 21, 20 St, Kaveh Industrial City, Saveh, Iran

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